

MELA CONTROL

BRIGHTENING • CORRECTING • PIGMENTATION CONTROL

A corrective bio-active serum that evens skin tone and corrects dark spots

MELA CONTROL delivers a brightening effect to the skin preventing excessive pigmentation. Fading the appearance of age spots, melasma, acne and uneven skin tone.

Reducing visible age spots, freckles and sun damage, along with the formation of post-inflammatory hyperpigmentation from acne.

It contains high concentration of **Niacinamide**, **Ascorbic acid**, **Tranexamic Acid**, and **Glutathione** to inhibit the transfer of melanin in skin cells, leading to a more even skin tone.



HOW TO USE

Like all other RETRO EV bioscience formulas, MELA CONTROL is designed to be combined with non-invasive clinical procedures such as microneedling or electroportation therapy. It can also be used as a topical serum post-laser.

1	Withdraw 2.5ml of MELA CONTROL using a syringe or attach the dropper cap included in the box
2	Apply the MELA CONTROL solution to one treatment zone at a time before going over with your treatment tool of choice, ensuring to massage manually with fingertips for absorption
3	Allow to dry before applying your chosen post-treatment cream

To maximise results, we recommend patients complement their in-clinic treatments with at-home routines.